OSTRZEŻENIA!

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WARNINGS!

1. The product is intended for children over the age of 8. The kit components have sharp edges and are made of materials that may be hazardous to children's safety.
2. To be used under the direct supervision of an adult.
3. Please keep the package and/or instruction manual. They contain important information that may be useful in the future.
4. Do not use this product for any purpose other than that for which it was designed.
5. The battery should not be recharged if it is not designed for this purpose.
6. Insert new batteries in the designated places, making sure that the battery terminals (+/-) are facing the right way.
7. Only one type of battery should be used at a time. Do not insert different types of batteries.
8. Batteries should be correctly seated in the battery compartment.
9. Remove depleted batteries immediately.
10. Do not use old batteries together with the new ones.
11. Never short-circuit the battery terminals. Do not short-circuit the battery connector.
12. Follow the disposal instructions on the battery.
13. Improper use of batteries can cause damage to them and may corrode the area around the battery posing a risk of fire, explosion, or damage to the device or health.

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Instructions for the kit Science Suitcase. Teaching aids with worksheets and lesson plans. Magnetism

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■ Age: 8+

The contents of the case allow the pupils to conduct basic magnetism experiments. The students have the opportunity to build electromagnets and test the magnetism of various materials. Thanks to the kit, the teacher can also demonstrate how magnets and magnetic field work. It is also possible to experience magnetic levitation. The lesson plans include the objectives, methods, forms and means necessary to introduce interesting learning activities. Additionally, they contain detailed descriptions of the lessons. The worksheets make it easier to carry out the experiments thanks to precise instructions.

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
Contents of the kit


- suitcase
dim. 55 x 15.5 x 46 cm
with accessories:

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Experiment #11

1. The students prepare the following materials:
 - electromagnetic device
 - wires
 - battery compartment
 - ferrous objects
 - objects that are not attracted to a magnet
2. Magnetism is the property that magnets have. It is a complex and uncontrolled process. This experiment allows the students to see if there is another way of creating magnetic materials.

Students install electromagnets. They then use the installed electromagnet to attract the ferromagnetic material.

What we can observe: a powered coil with an iron core attracts ferrous material. Conclusion: the energised coil is magnetic and has the same properties as a magnet. This device is called an electromagnet.

Now, the pupils construct a makeshift electromagnetic crane using an electromagnet to see which materials will be attracted to it.

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Cleaning and storing models

1. Keep away from contact with sharp, hard objects that could cause scratches.
2. Clean with a soft slightly damp cloth. Do not use sharp brushes, sponges, or scouring pads.
3. Avoid prolonged exposure to sunlight.
4. Use the kit's elements in a location where they will not be subjected to a fall or impact.
5. Store all kit components in the original suitcase in appropriate compartments.

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Experiment #8

1. The students prepare the following materials:

- bar magnet
- U-shaped magnet
- ferrous objects
- objects that are not attracted to a magnet

2. An experiment called “Magnetic Pole of a Magnet.”

We will use bar magnets to attract paper clips or pins. The two ends of a magnet are called magnetic poles. Now, we will do the same thing but with a large U-shaped magnet. This time, the two ends will also easily attract paper clips or pins. These two ends are the magnetic poles of the U-shaped magnet.

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Experiment #9

1. The students prepare the following materials:

- column magnet
- compass

2. Let's now learn how to use small magnetic needles. A small magnetic needle is the equivalent of a bar magnet due to the different directions that it is pointing. The north end is called the “north pole” (also called the “N” pole – north) and the leading end is called the south pole (also called the “S” pole – south).

Now use a small magnetic needle to determine the magnetic polarity of a magnet. If a cylindrical or bar magnet attracts the north pole of the small needle, this end is the south pole of that magnet. By doing so, we can determine the magnetic properties of the two magnetic poles of the magnet – North and South.

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Experiment #10


1. The students prepare the following materials:

- small magnetic balls
- magnetic hand horseshoe


2. The “Maglev” experiment. Stacking magnetic balls on a pole. Students will learn how objects can be suspended in the air.

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3 lesson plans in English and Polish on the following topics:

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24 single-sided worksheets, dim. 11 x 11 cm

The Magneticity Science Suitcase Worksheets are perfect for teaching about magnetism, magnets, electromagnets, and more. This comprehensive set of step-by-step experiment pictures contains all essential details needed by the teacher to help the student develop their scientific competence.

It consists of 11 experiments. Each labeled with a different colour. Every experiment is preceded by an image showing necessary materials. The next steps of the experiment are numbered. The pupils use the worksheets to improve their knowledge and practice their skills.

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Experiment #1

1. The students prepare the following materials:
 - ferrous materials
 - bar magnet
 - U-shaped magnet
2. The experiment is designed to show that magnets attract a large number of iron parts, which indicates how strong the magnet is. The strongest magnetic part of a magnet is called the magnetic pole.

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Experiment #2

1. Students prepare the following materials:
 - bar magnet
2. The experiment is designed to demonstrate that magnets have different magnetic poles and which we refer to as „N” and „S” poles respectively. Magnetic poles with the same name repel each other and magnetic poles with different names attract each other.

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Experiment #3

1. Students prepare the following materials:
 - magnetic bracket
 - ring magnets
2. The experiment is designed to show that, thanks to the principle of mutual repulsion of magnetic poles of the same name, it is possible to suspend a magnetic ring in the air.

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
Experiment #4

1. The students prepare the following materials:
 - carriages
 - ring magnets
2. The experiment shows that, because of the fact that magnetic poles with the same name repel each other and magnetic poles with different names attract each other, it is possible for trolleys to repel each other or for one trolley to pull the other.

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
Experiment #5

1. Students prepare the following materials:
 - electromagnetic device
 - wires
 - battery compartment
 - ferrous objects
2. The experiment shows that by connecting the device mentioned above to a battery, the coil is magnetised and can attract iron. This device is called an electromagnet.

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Experiment #6

1. The students prepare the following materials:
 - plastic labyrinth board
 - hand-held magnet
 - magnetic figures
2. An experiment called “Plastic Dolls in a Maze.” The pupils move plastic figures (with magnets) across the plastic board. The teacher asks them: Is this force a magnetic force?

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Experiment #7

1. The students prepare the following materials:
 - bar magnet
 - U-shaped magnet
 - iron filings
2. Students use magnets to attract iron filings. When the magnet gets closer to the iron filings, the black magnetic filings rise vertically and follow the movement of the magnet.